

The OPAL 480 s a fully hydraulic, water cooled hot mounting press for mounting samples before polishing. The process sequence is executed automatically. The fully automatic sliding closure allows quick and easy opening/closing of the mould assembly. Mould assemblies can be exchanged quickly and without tools. Also rectangular moulds can be produced. The robust aluminium construction is powder coated. With its solid construction and ATM hydraulic system, the operating noise of the machine is minimized.







SLIDING CLOSURE

The new sliding closure is a unique innovation and allows quick and easy opening/closing of the mould assembly. By actuating the lever, the closure slides to the correct position. After initiating the cycle the mould closes and reopens automatically on completion.

SIMPLE OPERATION

After pressing the start key, the fully automatic program sequence is executed with the preselected parameters.

The parameters can be changed and adjusted during the mounting process.



MOULD ASSEMBLIES

Our comprehensive selection of mould assemblies from \varnothing 50 - 70 mm resp. 30x60/40x60 mm (rectangular form) enables optimum monting of your sample. Two mounts can be produced simultaneously with the use of a spacer.

ELECTRONIC CONTROL

With its clear LCD operation panel you can easily control the hot mounting press. You can load and save eighteen different mounting processes.

HEATING

Our high-power heating system gives a very short warm-up time.
For delicate samples the machine can

For delicate samples the machine cardiust the heating curve infinitely.



COOLING

Not only the cooling time, but also the resulting cooling curve can be adjusted infinitely according to the mounting requirements. Thereby, a considerably lower water consumption is achieved for longer cooling times.

PRESSURE MODE

Besides the standard operation, you can choose between three pressure methods.

- » Pressure increases parallel with the pre-heating phase (ideal for hollow specimens, e.g. rings).
- » Pressure increases after the preset temperature is reached to produce mounts as gap free as possible (ideal for complex geometries).
- » When mounting highly filigree and pressure sensitive parts, the increase in pressure against the rising temperature is a big advantage.







BASIC MODULE



OPAL 480

Order No.: M0790000

- » Automatic hot mounting press
- » Fully hydraulic operation with automatic water cooling (with water saving function)
- » Completely new and simple automatic closure system
- » Computer controlled with automatic cycle processing
- » Four different pressure modes
- » Aluminium case, powder coated
- » Easily interchangeable moulds
- » Rectangular moulds

EQUIPMENT

POWER SUPPLY (select equipment)

Equipment 1

230 V/50 Hz (1/N/PE) Order No: A0790000 Equipment 2

110 V/60 Hz (1/N/PE) Order No.: A0790001

ACCESSORIES

COOLING UNIT

Recirculation unit



» Cooling» Motor» Flow pressure» Tank0.81 kW0.6 kW2.8 bar6 Ltr.

» WxHxD 440x345x565mm

System cooling unit

» Integrated in drawer

» Cooling» Motor» Flow pressure» Tank0.81 kW0.6 kW2.8 bar6 Ltr.

230 V/50 Hz (1/N/PE) Order No.: Z0730007 110 V/60 Hz (1/N/PE)

Order No.: Z0730008

230 V/50 Hz (1N/PE) Order No.: A5800118

110 V/60 Hz (1/N/PE) Order No.: A5800141



Water cooling with water saving-function

ACCESSORIES

MOULD ASSEMBLIES (CYLINDER, UPPER- AND LOWER RAM)

	Order No.:
Ø 50 mm	Z0790000
Ø 60 mm	Z0790001
Ø 70 mm	Z0790002
\square 30 x 60 mm	Z0790003
□ 40 x 60 mm	Z0790004

SPACER FOR DOUBLE MOUNTS

Order No.:
Ø 50 mm
00702118
Ø 60 mm
00702025
Ø 70 mm
00702026



CONNECTION SET

Order No.: Z0730003 » 2 pressure tubes R ½"; 2 m lang

» Adapter for ¾"

Order No.: Z0730004

» Connection set for two Opal presses




SYSTEM LAB CABINET


Order No.: M5800048

» WxHxD 880x800x800 mm

» System lab assembly

SPECIFICATIONS





Moulds	Ø 50 - 70 mm
Moulds (rectangular)	30x60/40x60 mm
Connected load	3.6 kVA
Power supply	1 230 V/50 Hz (1/N/PE)
	2 110 V/60 Hz (1/N/PE)
Heating	4 x 630 Watt
Temperature range	0 – 200 °C
Heating time	variable / 0 - 30 min
Cooling time	variable / 0 - 30 min
Pressure	up to 300 bar depending on mould
WxHxD	390 x 450 x 500 mm
Weight	68 kg