



Combination for 2-lines

SYSTEM-INTEGRATION

The Opal 460 and 480 are available not only as single units but also as a bench standing unit. By integrating the mounting presses into our System, we produce a comfortable working height.

Optionally, the SYSTEM-MOUNTING PRESS can be equipped with a closed cooling water system, so it becomes independent

from any mains water supply. In addition, the water consumption is reduced to a single filling (6 litres).

The system-mounting presses are available as twin- or triple powder-coated aluminum units and fit perfectly into our ATM System Lab.



SYSTEM-HOT MOUNTING PRESS



- » automatic hot mounting press in System Lab design, optionally available with Opal 460, 480 and HPL-tabletop
- » every mounting press operates independently
- » automatic water cooling (with water-miser function)
- » automatic sliding closure system
- » fully hydraulic operation
- » computer controlled with automatic cycle system
- » three different pressure modes
- » easy change of mould assemblies
- » made of powder-coated aluminium
- » e-Lab ready

| | |
|---------------|---------------|
| 2-lines WxHxD | 88x110x80 cm |
| 3-lines WxHxD | 132x110x80 cm |