

FIXATION SYSTEM FOR SILICON GRINDING PAPERS WITHOUT ADHESIVE/FOIL BACKING

Qprep Clamping Ring

The QPREP clamping ring is made of plastic and allows grinding papers without adhesive / foil backing to be attached directly to the working disc. This ensures very low resilience and minimizes edge rounding during the grinding process. For changing the grinding paper, the clamping ring must be re-attached each time.



Item No.	Unit	Description
CLAMPING RING FOR GRINDING PAPERS		
Ø		
• made from plastic		
05610228	1 Pc.	200 mm
05610227	1 Pc.	230 mm
05610226	1 Pc.	250 mm
05610225	1 Pc.	300 mm

FIXATION SYSTEM FOR SILICON GRINDING PAPERS WITHOUT ADHESIVE/FOIL BACKING

Qprep Quick-Tap

The reusable fast and time saving fixation system for plain backed grinding paper. As an alternative to the clamping ring and adhesive film, the innovative QPREP Quick-Tap is ideally suited for fixing grinding paper without adhesive/foil backing. Designed to be reusable, the thin layer of our new innovative fixation system resists high shear forces and ensures that the plain backed paper adheres securely to the Quick-Tap. Like all new Galaxy products, the Quick-Tap features the Qprep Antislip backside for secure retention of the Qprep magnetic film.


PRODUCT ADVANTAGES

- | Optimized adhesion
- | Low impact elasticity
- | Reduction of pencil effect and edge rounding
- | Possibility to position the specimen holder beyond the edge of the grinding paper during semi-automatic and automatic specimen preparation
- | Sustainable and resource-saving compared to using grinding paper with adhesive/foil backing
- | Allows hundreds of paper changes

RECOMMENDED APPLICATIONS

- | For use with grinding paper without adhesive/foil backing
- | Quick and easy as well as repeated fixation of grinding paper
- | For direct attachment to magnetic foil

Item No.	Unit	Description
QUICK-TAP		
Ø		
GALAXY Quick-Tap (for using with magnetic foil)		
95017484	1 Pc.	250 mm
95017485	1 Pc.	300 mm